**NOTICE D'UTILISATION APPLICATEUR SPIRE'IT®**  
**DIRECTIONS FOR USE SPIRE'IT®APPLIER**  
**INSTRUÇÕES DE USO APLICADOR SPIRE'IT®**  
**GEBRAUCHSANWEISUNG APPLIKATOR SPIRE'IT®**  
**ISTRUZIONI PER L'USO APPLICATORE SPIRE'IT®**  
**WSKAZÓWKI DOTYCZĄCE STOSOWANIA APLIKATORA SPIRE'IT®**

**À LIRE ATTENTIVEMENT** par l'opérateur ainsi que par le personnel de bloc opératoire  
**TO BE CAREFULLY READ AND STRICTLY OBSERVED** by the surgeon and the staff of the operating theatre.  
**LEA DETENIDAMENTE** por el operador así como por el personal del bloque quirúrgico.

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This instruction for use is not a reference document on surgical stappling or suturing but an information about the use of SpiRE'IT® suture applicator.  
 The use of this device has to be carried out exclusively in operating rooms with requested aseptic conditions mentioned for a surgical procedure.  
 The use of this device is strictly reserved to specialized surgeons, sensitized to the use of SpiRE'IT®. A video is available on [www.microval.fr](http://www.microval.fr).  
 The non-observance of these recommendations may contribute to a reduced life duration of the device. The surgeon is responsible for any complication that may result from erroneous indication, incorrect operative technical or inadequate aseptic precautions.  
 The patient must be informed of the precautions that he will have to take after the implantation as well as the contra-indications and undesirable side-effects that may occur.  
 Clinical results are depending on the selected material, the surgical indication and the quality of the surgical operation.

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**INDICATIONS :**  
 MICROVAL's SpiRE'IT® suture applicator is used in general endoscopic or open surgery for meshes fixation, during hernia repair and approximation of tissues in various surgical specialties.  
**DESCRIPTION OF THE DEVICE :**  
 MICROVAL's SpiRE'IT® suture applicator is a sterile device designed for a multiple use during one same operation.  
 Applicators for endoscopy are to be used with appropriate sized trocar sleeve diameter 5mm minimum.  
 MICROVAL's SpiRE'IT® suture applicer allows to implant memory shape nickel-Titanium alloys (NITINOL) having a double spire 4mm diameter shape.  
 MICROVAL's SpiRE'IT® suture applicer is delivered sterile ready to use, with 10 suture spires.  
 The applicer implant one fixation at each pressure on the trigger (A). At expulsion of its place, the spire closes itself in spire suture shape, fixes or approximates tissues and/or meshes.  
 (A) Trigger, (B) extremity, (C) cartridge inside B, (D) rotation indicator at the rear of the applicer.

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
**INSTRUCTIONS FOR USE :**  
 Check good compatibility of all instruments to each other.  
 Remove MICROVAL's SpiRE'IT® suture applicer of its packaging according usual techniques for sterile products.  
 To avoid any damage, do not drop the applicer on the sterile area.  
 During endoscopic operation, introduce the applicer through an adapted sleeve trocar (5mm diameter minimum).  
 Place the tip of the applicer at the center of the operating site and perpendicular to tissues or mesh.  
 In the case of approximation of tissues, use a grip forceps before any spire implementation.  
 Apply light pressure to the tissue site during SpiRE'IT® suture placement.  
 During placement of a spire internally, apply an external counter-pressure to the area being fixed when possible to ensure secure suture placement.  
 To fire the applicer, simultaneously press the tip against the tissue or mesh and squeeze completely the trigger (A) and release to set the instrument for the next use.  
 Do not release the trigger roughly during applying a suture.  
 Do not apply excessive force when squeezing the trigger the applicer may jam.  
 If a SpiRE'IT® suture becomes jammed remove all the applicer.  
 Check the fixed area to ensure the good shape and place of the SpiRE'IT®.  
 Repeat all the above steps along the length of the area to be approximated and/or fixed until the procedure is complete.  
**SINGLE USE CARTRIDGE (C) :**  
 A rotation indicator (D) located at the rear of the applicer allows to spot on the condition of the cartridge. Check the initial position of this indicator at the beginning of the procedure. Once this indicator made one 360° rotation, the 10 fixation spires have been implanted.

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**PRECAUTIONS OF USE :**  
 An excessive pressure on the tip of the applicer during firing may jam and block it.  
 After implantation of any SpiRE'IT® suture with the applicer, always ensure carefully the haemostasis. It is then possible to care the bleeding with electro-coagulation or manual suture.

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**DISCARD AFTER USE**  
 All contaminated applicer has to be decontaminated if it has to be return to MICROVAL.

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**CONTRA-INDICATIONS :**  
 Do not use SpiRE'IT® suture applicer on vascular, nervous tissues nor on hard bodies (bones).  
 It is however possible to place SpiRE'IT® sutures above underlying bones, vessels or viscera, but there must be a distance of not less 4.5 mm from the surface of the tissue to the underlying structure.  
 Do not use SpiRE'IT® suture applicer to perform haemostasis.  
 Do not use SpiRE'IT® suture applicer on a site where it is impossible to visually check haemostasis.  
 Store at room temperature  
 Avoid prolonged exposure to elevated temperature.

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**ADVERSE REACTIONS:**  
 Seroma/hematoma formation; Urinary retention; Neuralgias; Testicular pain and swelling; Mesh infection and wound infection; Recurrence.

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